



WhaleTeq MECG 2.0 SDK

Version: 1.0.3.0 Build: 20201104

=====

1. redesign periodical waveform (sine/ triangle/ pulse); load as continuous waveform instead of 10-seconds waveform

Version: 1.0.2.0 Build: 20201015

=====

1. update ECG_HEADER structure; field ECG_SIGNAL Signal[1] is defined as a variable-length array
2. add MEGConnect() to connect to the device

Version: 1.0.1.2 Build: 20180801

=====

1. added support new MEGC HW (DAC 16bit)
2. added support to output sine/triangle wave with different amplitudes to different leads

Version: 1.0.0.9 Build: 20170925

=====

1. initial version