



WhaleTeq SEEG 100 PC Software

Version: 1.0.5.10 Build: 20190906

=====

Fix bugs in SEEG-100E Assistant 80601-2-26:2019.

Version: 1.0.5.9 Build: 20190802

=====

Support SEEG-100E(Assistant 80601-2-26:2019).

Version: 1.0.5.8 Build: 20190313

=====

Fix first time database downloading failed bug.
Change limit of "Variable DC Offset" in Calibration from 5% to 1%.
Fix "Live Firmware Update" notification failed bug.
Fix Serial Number display error on windows bar.

Version: 1.0.5.6 Build: 20180309

=====

Live Update for Software and Firmware.
Change the initial value of amplitude and frequency to "100uV" and "10Hz".
Refine "Load EDF Database" and "Raw Data Replay".
Installation version with Certification supported.
Support "Epileptic Database" - [chbmit\(https://physionet.org/physiobank/database/chbmit/\)](https://physionet.org/physiobank/database/chbmit/)

Version: 1.0.0.3 Build: 20170209

=====

Add ISP (In System Programming) Feature, it includes below two functions:
- Firmware Update
- Binding Assistant active key to SEEG device.
The feature is only supported with new hardware and firmware.



RELEASE NOTE

TEL | +886-2-2596-0701
FAX | +886-2-2596-0702
eMail | service@whaleteq.com

Version: 1.0.0.2 Build: 20160519

=====

New Feature: Support EDF file and File Manager

Version: 1.0.0.1 Build: 20151028

=====

Initial release version for EEG AP.